



DETECT & DISSOLVE
BIOFILM MATRIX

Instructions for use

Use enziDent® each time instruments are cleaned for complete removal of soil and to prevent build-up of incrustated soil and biofilm matrix.

- A.** Standard dosage : **1%***. Increase dosage to **2%*** in case of hard water or heavy soiling or temperature of water < 30°C
- B.** Contact time : **15 minutes**
- C.** Ideal temperature : **40 to 45°C** (minimum 30°C and maximum 55°C)

*: See dilution table below

More details

- Instruments should be rinsed prior to immersion in enziDent® only if product used for initial soak has an extreme pH (>10) or contains oxidizing agents (peracetic acid, hydrogen peroxide or chlorine-based products) or phenols.
- Respect recommended temperatures for optimal performance; efficacy is not guaranteed > 55°C.
- Rinse abundantly with water before disinfection and sterilization.
- Enzymatic activity is maintained for 8h following dilution.
- Renew baths frequently according to soil levels.

Delivery units

Product reference	Packaging
OL20806	5L
OL20807	2 x 1L
OL20808	6 x 1L
OL20809	2 x 5L

Dilution table (Dosing head: 25ml)

Bath volume (l)	0,5	1	2	3	4	5	10	15
enziDent® 1%-2% (in ml)	5 - 10	10 - 20	20 - 40	30 - 60	40 - 80	50 - 100	100 - 200	150 - 300

Chemical Properties

Appearance :	Yellow liquid
pH concentrate :	8.5
pH solution 1%:	7.8 +/-0.5
Density:	1.053+/-0.01

The color is natural and may vary from one batch to another with no impact on the performance of the product

Main components

1 - 5% non-ionic surfactants;
 1 - 5% anionic surfactants;
 <1% sequestrants;
 <5% corrosion inhibitors;
 Patented multi-enzymatic compound

Precautions

- **Wear gloves.** In case of prolonged use, gloves covering forearms are recommended. Refer to full safety data sheet before use.
- **Store in the original container, closed, between +4°C and +25°C.** For optimal performance respect the « Best Before » date on the label.
- **End of life:** an agreed collector must eliminate unused product.

Maahantuonti & markkinointi:

SYMPPIIS International Oy
 PL 67
 02771 Espoo
 020 735 3333 - info@symppis.com

www.symppis.com